




## Product specification

**Tip comb, free from DNA/DNase/RNase/PCR inhibitors, PP**



## Product information

Product number	82.1972.800
----------------	-------------

## Product descriptions

Name	Tip comb, free from DNA/DNase/RNase/PCR inhibitors, PP
Product number	82.1972.800
Product description	Tip comb, compatible with KingFisher™ Flex/Presto, Bio Sprint 96, free from DNA/DNase/RNase/PCR inhibitors, material: PP, 2 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Diameter	7.52 mm
Width of product	85 mm
Length of product	127.76 mm
Height of product	41.85 mm

## Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	no medical device   IVD
Purity standard	free from DNA/DNase/RNase/PCR inhibitors
Batched	yes

## Packaging

Minimum order quantity	20
Piece(s) / case	20
Piece(s) / pallet	2160
Height of case	195 mm
Width of case	223 mm
Depth of case	308 mm
Volume of case	0.0134 cbm
Weight of case	1.2776 kg
EAN inner box	4038917958080
EAN case	4038917958097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT