




## Product specification

### Tip comb, PCR Performance Tested, PP



## Product information

Product number	82.1972.800
----------------	-------------

## Product descriptions

Name	Tip comb, PCR Performance Tested, PP
Product number	82.1972.800
Product description	Tip comb, compatible with KingFisher™ Flex/Presto, Bio Sprint 96, PCR Performance Tested, material: PP, 2 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	7.52 mm
Width of product	85 mm
Length of product	127.76 mm
Height of product	41.85 mm

## Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Purity standard	PCR Performance Tested
-----------------	------------------------

Batched	yes
---------	-----

## Packaging

Piece(s) / case	20
-----------------	----

Piece(s) / pallet	2160
-------------------	------

Height of case	195 mm
----------------	--------

Width of case	223 mm
---------------	--------

Depth of case	308 mm
---------------	--------

Volume of case	0.0134 cbm
----------------	------------

Weight of case	1.4 kg
----------------	--------

EAN inner box	4038917958080
---------------	---------------

EAN case	4038917958097
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT