



Product specification

Micro sample tube Citrate 9NC, 1 ml, push cap, ISO



Product information

Product number	41.1506.100
----------------	-------------

Product descriptions

Name	Micro sample tube Citrate 9NC, 1 ml, push cap, ISO
Product number	41.1506.100
Product description	Micro sample tube Citrate 9NC, venous blood collection, preparation: Citrate 9NC, nominal volume: 1 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: blue, colour code ISO, with plastic label, label/print: blue, skirted conical base, 100 piece(s)/polystyrene container

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	blue

Size

Diameter	10.8 mm
Length excluding cap	40 mm
Length including cap	44 mm
Length of product	44 mm

Sample volume	1 ml
---------------	------

Product characteristics

Base shape	skirted conical base
Label/Print	With plastic label
Color of print/label	blue
Graduation	yes
Preparation concentration	0.106 mol/l
Cap	flat cap
Closure type	Push cap
Type of collection	venous
Centrifugation recommended	1500 x g
Centrifugation time	10 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	44000
Height of case	155 mm
Width of case	348 mm
Depth of case	568 mm
Volume of case	0.0306 cbm
Weight of case	2.44 kg
EAN inner box	4038917566032
EAN case	4038917566025

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT