




Product specification

Pipette tip, 10 ml, blue, 100 piece(s)/bag



Product information

Product number	70.1187.002
----------------	-------------

Product descriptions

Name	Pipette tip, 10 ml, blue, 100 piece(s)/bag
Product number	70.1187.002
Product description	Pipette tip, working volume: 10 ml, blue, filling level rings, suitable for Sarpette® M, 100 piece(s)/bag

Material & colours

Color of product	blue
------------------	------

Size

Length of product	144 mm
Max. volume	10 ml

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	20000
Height of inner box	65 mm
Width of inner box	325 mm
Depth of inner box	520 mm
Height of case	269 mm
Width of case	381 mm
Depth of case	457 mm
Volume of case	0.0468 cbm
Weight of case	5.4 kg
EAN inner box	4038917951531
EAN case	4038917951548

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT