




## Product specification

**Pipette tip, 10 ml, blue, 25 piece(s)/box**



### Product information

Product number	70.1187.102
----------------	-------------

### Product descriptions

Name	Pipette tip, 10 ml, blue, 25 piece(s)/box
Product number	70.1187.102
Product description	Pipette tip, working volume: 10 ml, blue, filling level rings, suitable for Sarpette® M, 25 piece(s)/box

### Material & colours

Color of product	blue
------------------	------

### Size

Length of product	144 mm
Max. volume	10 ml

### Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

### Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	25
Piece(s) / case	250
Piece(s) / pallet	7500
Height of inner box	148 mm
Width of inner box	117 mm
Depth of inner box	150 mm
Height of case	152 mm
Width of case	303 mm
Depth of case	603 mm
Volume of case	0.0278 cbm
Weight of case	3.6 kg
EAN inner box	4038917951494
EAN case	4038917951500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT