



Product specification

Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugati



Product information

Product number	53.419
----------------	--------

Product descriptions

Name	Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation
------	--

Product number	53.419
----------------	--------

Product description	Seraplas® valve filter, V13, blue, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case
---------------------	--

Material & colours

Color of product	blue
------------------	------

Size

Length of product	92 mm
-------------------	-------

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched	no
---------	----

Safety data sheet	no
-------------------	----

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	75000
Height of inner box	98 mm
Width of inner box	99 mm
Depth of inner box	119 mm
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm
Weight of case	2.15 kg
EAN inner box	4038917081481
EAN case	4038917000574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT