



Product specification

Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugati



Product information

| | |
|----------------|--------|
| Product number | 53.419 |
|----------------|--------|

Product descriptions

| | |
|---------------------|--|
| Name | Seraplas® valve filter, blue, for separation of serum/plasma from the blood cells after centrifugation |
| Product number | 53.419 |
| Product description | Seraplas® valve filter, V13, blue, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case |

Material & colours

| | |
|------------------|------|
| Color of product | blue |
|------------------|------|

Size

| | |
|-------------------|-------|
| Length of product | 92 mm |
|-------------------|-------|

Purity & certification

| | |
|-------------------|---------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |
| Safety data sheet | no |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 100 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 75000 |
| Height of inner box | 98 mm |
| Width of inner box | 99 mm |
| Depth of inner box | 119 mm |
| Height of case | 140 mm |
| Width of case | 293 mm |
| Depth of case | 513 mm |
| Volume of case | 0.021 cbm |
| Weight of case | 2.15 kg |
| EAN inner box | 4038917081481 |
| EAN case | 4038917000574 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT