



## Product specification

**Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP**



### Product information

Product number	60.550.020
----------------	------------

### Product descriptions

Name	Screw cap tube, 7 ml, (LxØ): 82 x 13 mm, PP
Product number	60.550.020
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	13 mm
Length of product	82 mm
Height of product	82 mm

### Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Screw cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	4000 x g

## Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	3000
Piece(s) / pallet	48000
Volume of case	0.0813 cbm
Weight of case	10.251 kg
EAN inner box	4038917112284
EAN case	4038917032537

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT