




Product specification

**Transfer pipette, 1 ml, (LxW): 104 x 10 mm, LD-PE,
transparent**



Product information

Product number	86.1180
----------------	---------

Product descriptions

Name	Transfer pipette, 1 ml, (LxW): 104 x 10 mm, LD-PE, transparent
Product number	86.1180
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 104 x 10 mm, material: LD-PE, transparent, droplet size: 17 µl - 22 µl, 750 piece(s)/case

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Color of product	transparent

Size

Width of product	10 mm
Length of product	104 mm
Sample volume	1 ml
Droplet size min-max	17 µl - 22 µl

Product characteristics

Graduation	no
------------	----

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	3000
Piece(s) / inner box	750
Piece(s) / case	3000
Piece(s) / pallet	252000
Height of inner box	121 mm
Width of inner box	126 mm
Depth of inner box	261 mm
Height of case	140 mm
Width of case	263 mm
Depth of case	533 mm
Volume of case	0.0196 cbm
Weight of case	2.07 kg
EAN inner box	4038917085649
EAN case	4038917004282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT