



Product specification

Screw cap, natural, suitable for 25ml , 30ml tubes



Product information

Product number	65.9923.335
----------------	-------------

Product descriptions

Name	Screw cap, natural, suitable for 25ml , 30ml Tubes
Product number	65.9923.335
Product description	Screw cap, natural, suitable for 25ml , 30ml tubes, 250 piece(s)/bag

Material & colours

Color of product	natural
------------------	---------

Purity & certification

Pyrogen-free	no
Latex-free	yes
Safety data sheet	no

Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	40000
Height of case	158 mm
Width of case	224 mm

Depth of case	271 mm
Volume of case	0.0099 cbm
Weight of case	1.4 kg
EAN inner box	4038917153874
EAN case	4038917079013

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT