



Product specification

Sample tube, suitable for Boehringer ES 300® and ES 600®, transparent



Product information

Product number	73.663
----------------	--------

Product descriptions

Name	Sample tube, suitable for Boehringer ES 300® and ES 600®, transparent
Product number	73.663
Product description	Sample tube, suitable for Boehringer ES 300® and ES 600®, material: PS, 2 ml, (LxØ): 40 x 12 mm, transparent, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	12 mm
Length of product	40 mm
Height of product	40 mm
Diameter of opening	10 mm
Max. volume	2 ml

Product characteristics

Base shape	Conical base
------------	--------------

Graduation	no
------------	----

Closure type	Push cap
--------------	----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE - manufacturer's self-declaration
--------------	--------------------------------------

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	500
----------------------	-----

Piece(s) / case	5000
-----------------	------

Piece(s) / pallet	180000
-------------------	--------

Height of case	310 mm
----------------	--------

Width of case	398 mm
---------------	--------

Depth of case	398 mm
---------------	--------

Volume of case	0.0491 cbm
----------------	------------

Weight of case	6.72 kg
----------------	---------

EAN inner box	4038917083195
---------------	---------------

EAN case	4038917002158
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT