



Product specification

Sample tube, suitable for Abbott analyser AXSYM, transparent



Product information

| | |
|----------------|--------|
| Product number | 73.667 |
|----------------|--------|

Product descriptions

| | |
|---------------------|---|
| Name | Sample tube, suitable for Abbott analyser AXSYM, transparent |
| Product number | 73.667 |
| Product description | Sample tube, suitable for Abbott analyser AXSYM, material: PS, 2.5 ml, (LxØ): 53 x 12 mm, transparent, 500 piece(s)/bag |

Material & colours

| | |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Color of product | transparent |

Size

| | |
|---------------------|--------|
| Diameter | 12 mm |
| Length of product | 53 mm |
| Height of product | 53 mm |
| Diameter of opening | 10 mm |
| Max. volume | 2.5 ml |

Product characteristics

| | |
|------------|--------------|
| Base shape | Conical base |
| Graduation | no |

Purity & certification

| | |
|------------------|-----------------------|
| Product category | diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 2500 |
| Piece(s) / inner box | 500 |
| Piece(s) / case | 2500 |
| Piece(s) / pallet | 90000 |
| Height of case | 265 mm |
| Width of case | 398 mm |
| Depth of case | 398 mm |
| Volume of case | 0.042 cbm |
| Weight of case | 6.78 kg |
| EAN inner box | 4038917083232 |
| EAN case | 4038917002189 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT