



Product specification

Sample tube, suitable for Abbott analyser AXSYM, transparent



Product information

Product number	73.667
----------------	--------

Product descriptions

Name	Sample tube, suitable for Abbott analyser AXSYM, transparent
Product number	73.667
Product description	Sample tube, suitable for Abbott analyser AXSYM, material: PS, 2.5 ml, (LxØ): 53 x 12 mm, transparent, 500 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	12 mm
Length of product	53 mm
Height of product	53 mm
Diameter of opening	10 mm
Max. volume	2.5 ml

Product characteristics

Base shape	Conical base
------------	--------------

Graduation

no

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2500
Piece(s) / pallet	90000
Height of case	265 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.042 cbm
Weight of case	6.78 kg
EAN inner box	4038917083232
EAN case	4038917002189

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT