



Product specification

Snap-on lid, with spout, PE, red, 100 piece(s)/bag



Product information

Product number	76.1477.100
----------------	-------------

Product descriptions

Name	Snap-on lid, with spout, PE, red, 100 piece(s)/bag
Product number	76.1477.100
Product description	Snap-on lid, with spout, suitable for 125-ml container 75.1355, material: PE, red, 100 piece(s)/bag

Material & colours

Product material	Polyethylene (PE)
Color of product	red

Product characteristics

Design	with spout
Closure type	Snap-on lid

Purity & certification

Product category	diagnostic device, CE
CE certified	CE supplier
Batched	yes

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	24000
Height of case	590 mm
Width of case	265 mm
Depth of case	390 mm
Volume of case	0.061 cbm
Weight of case	4.57 kg
EAN inner box	4038917129411
EAN case	4042836001604

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT