



## Product specification

**Absorbent liner, suitable for Mailing container 114 x 44 mm, (LxW): 135 x 95 mm**



## Product information

Product number	78.669
----------------	--------

## Product descriptions

Name	Absorbent liner, suitable for mailing container 114 x 44 mm, (LxW): 135 x 95 mm
Product number	78.669
Product description	Absorbent liner, suitable for mailing container 114 x 44 mm, (LxW): 135 x 95 mm, 100 piece(s)/case

## Material & colours

Product material	Cellulose
------------------	-----------

## Size

Width of product	95 mm
Length of product	135 mm

## Purity & certification

Product category	no medical device   IVD
Batched	no

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	56000
Height of case	160 mm
Width of case	385 mm
Depth of case	385 mm
Volume of case	0.0237 cbm
Weight of case	3.04 kg
EAN inner box	4038917083249
EAN case	4038917002196

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT