




Product specification

Screw cap, for tissue culture flask T-25, without filter



Product information

| | |
|----------------|-------------|
| Product number | 83.3990.025 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|--|
| Name | Screw cap, for Tissue culture flask T-25, without filter |
| Product number | 83.3990.025 |
| Product description | Screw cap, for tissue culture flask T-25, 2-position screw cap, without filter, red, TC Tested, 1 piece(s)/bag |

Material & colours

| | |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Color of product | red |
| Color of cap | red |

Product characteristics

| | |
|--------------|----------------------|
| Filter | no |
| Cap | 2-position screw cap |
| Closure type | Screw cap |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | TC Tested |
| Sterilisation | Gamma irradiation |

| | |
|--------------|-----|
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 100 |
| Piece(s) / inner box | 25 |
| Piece(s) / case | 100 |
| Piece(s) / pallet | 16800 |
| Height of case | 150 mm |
| Width of case | 218 mm |
| Depth of case | 288 mm |
| Volume of case | 0.0094 cbm |
| Weight of case | 0.57 kg |
| EAN inner box | 4038917359924 |
| EAN case | 4038917359917 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT