




## Product specification

### Screw cap, for tissue culture flask T-25, without filter



## Product information

Product number	83.3990.025
----------------	-------------

## Product descriptions

Name	Screw cap, for Tissue culture flask T-25, without filter
Product number	83.3990.025
Product description	Screw cap, for tissue culture flask T-25, 2-position screw cap, without filter, red, TC Tested, 1 piece(s)/bag

## Material & colours

Color of product	red
Color of cap	red

## Product characteristics

Filter	no
Cap	2-position screw cap
Closure type	Screw cap

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	TC Tested
Sterilisation	Gamma irradiation

Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	16800
Height of case	150 mm
Width of case	218 mm
Depth of case	288 mm
Volume of case	0.0094 cbm
Weight of case	0.57 kg
EAN inner box	4038917359924
EAN case	4038917359917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT