




## Product specification

# ESR rack, suitable for S-Sedivette®, 10 measuring stations



## Product information

Product number	90.1090
----------------	---------

## Product descriptions

Name	ESR rack, suitable for S-Sedivette®, 10 measuring stations
------	--

Product number	90.1090
----------------	---------

Product description	ESR rack, suitable for S-Sedivette®, 10 measuring stations
---------------------	--

## Product characteristics

Number of positions	10
---------------------	----

## Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE
--------------	----

Batched	yes
---------	-----

## Packaging

Piece(s) / case	16
-----------------	----

Piece(s) / pallet	320
-------------------	-----

Height of inner box	211 mm
---------------------	--------

Width of inner box	73 mm
Depth of inner box	73 mm
Height of case	350 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0833 cbm
Weight of case	11.75 kg
EAN inner box	4038917085304
EAN case	4038917004008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT