



Product specification

Tube segment opener, for bed side test cards, HD-PE, white



Product information

Product number	95.1000.100
----------------	-------------

Product descriptions

Name	Tube segment opener, for bed side test cards, HD-PE, white
Product number	95.1000.100
Product description	Tube segment opener, for bed side test cards, material: HD-PE, white, 200 piece(s)/case

Material & colours

Product material	High Density Polyethylene (HD-PE)
Color of product	white

Product characteristics

Design	for bed side test cards
--------	-------------------------

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	2000
Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	242000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	4.48 kg
EAN inner box	4038917124607
EAN case	4038917028585

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT