




Product specification

End-to-end capillary, K2 EDTA, 20 µl, Ø: 1.28 mm, glass



Product information

Product number	19.447.001
----------------	------------

Product descriptions

Name	End-to-end capillary, K2 EDTA, 20 µl, Ø: 1.28 mm, glass
Product number	19.447.001
Product description	End-to-end capillary, K2 EDTA, capillary blood collection, preparation: K2 EDTA, nominal volume: 20 µl, Ø: 1.28 mm, glass, 100 piece(s)/container

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	1.28 mm
Length of product	29 mm
Sample volume	20 µl
Blood volume	20 µl

Product characteristics

Type of collection	capillary
--------------------	-----------

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	50000
Piece(s) / pallet	2400000
Height of inner box	40 mm
Width of inner box	50 mm
Depth of inner box	117 mm
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	8.65 kg
EAN inner box	5706860174428
EAN case	5706860174435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT