



Product specification

Filter tip, 1,000 µl, black, 480 piece(s)/blister



Product information

Product number	70.1193.150
----------------	-------------

Product descriptions

Name	Filter tip, 1,000 µl, black, 480 piece(s)/blister
Product number	70.1193.150
Product description	Filter tip, working volume: 1,000 µl, black, suitable for application on SARSTEDT IVARO AP devices, 480 piece(s)/blister

Material & colours

Color of product	black
------------------	-------

Size

Length of product	95 mm
Max. volume	1000 µl

Product characteristics

Filter	yes
--------	-----

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

Packaging

Piece(s) / inner box	3840
Piece(s) / case	3840
Piece(s) / pallet	46080
Height of case	500 mm
Width of case	343 mm
Depth of case	563 mm
Volume of case	0.0966 cbm
Weight of case	7.8989 kg
EAN inner box	4038917434768
EAN case	4038917434751

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT