



Product specification

7 ml, (LxØ): 82 x 13 mm, PP



Product information

Product number	60.550.300
----------------	------------

Product descriptions

Name	7 ml, (LxØ): 82 x 13 mm, PP
Product number	60.550.300
Product description	working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, screw cap, without cap, sterile, 1,000 piece(s)/bag

Size

Diameter	13 mm
Length of product	82 mm
Height of product	82 mm

Product characteristics

Base shape	Flat base
Graduation	no
Closure type	Screw cap
Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	4000 x g

Purity & certification

Purity standard	sterile
-----------------	---------

Batched

yes

Packaging

Piece(s) / inner box	1000
Piece(s) / case	3000
Piece(s) / pallet	54000
Volume of case	0.0844 cbm
Weight of case	10.02 kg
EAN inner box	4038917112642
EAN case	4038917032551

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT