



Product specification

Screw cap, red, suitable for 15 ml tubes



Product information

Product number	65.797.001
----------------	------------

Product descriptions

Name	Screw cap, red, suitable for 15 ml tubes
Product number	65.797.001
Product description	Screw cap, red, suitable for 15 ml tubes, 1,000 piece(s)/bag

Material & colours

Color of product	red
Color of cap	red

Purity & certification

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	120000
Volume of case	0.0218 cbm
EAN inner box	4038917430470
EAN case	4038917430463

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT