



Product specification

Pipetting aids, Rotatable, storage capacity: 6, for Single-channel pipettes up to 1,000 µl, 1 piece(



Product information

Product number	95.1760.200
----------------	-------------

Product descriptions

Name	Pipetting aids, Rotatable, storage capacity: 6, for Single-channel pipettes up to 1,000 µl, 1 piece(s)/case
Product number	95.1760.200
Product description	Pipetting aids, Rotatable, storage capacity: 6, for Single-channel pipettes up to 1,000 µl, 1 piece(s)/case

Material & colours

Color of product	white/red
------------------	-----------

Product characteristics

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	144
-------------------	-----

Height of inner box	150 mm
Width of inner box	150 mm
Depth of inner box	230 mm
Height of case	260 mm
Width of case	337 mm
Depth of case	365 mm
Volume of case	0.032 cbm
Weight of case	3.158 kg
EAN inner box	4038917937863
EAN case	4038917937870

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT