



Product specification

PCR strip of 4, 200 µl, Biosphere® plus, transparent, PP, flat cap



Product information

Product number	72.990
----------------	--------

Product descriptions

Name	PCR strip of 4, 200 µl, Biosphere® plus, transparent, PP, flat cap
Product number	72.990
Product description	PCR strip of 4, up to 200 µl, High Profile, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 50 piece(s)/bag

Material & colours

Color of product	transparent
Profile	High Profile

Size

Diameter	5.35 mm
Width of product	10.6 mm
Length of product	37.2 mm
Height of product	20.65 mm
Max. volume	200 µl
Sample volume	0.2 ml

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Snap-on lid
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	200
Piece(s) / case	400
Piece(s) / pallet	57600
Height of inner box	97 mm
Width of inner box	151 mm
Depth of inner box	189 mm
Height of case	145 mm
Width of case	213 mm
Depth of case	323 mm
Volume of case	0.01 cbm
Weight of case	0.63 kg
EAN inner box	4038917084871
EAN case	4038917030045

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT