



## Product specification

### Container with screw cap, 100 ml, PP, transparent



## Product information

Product number	75.1354.001
----------------	-------------

## Product descriptions

Name	Container with screw cap, 100 ml, PP, transparent
Product number	75.1354.001
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø opening: 58 mm, transparent, graduated, material: PP, screw cap, cap: green, cap enclosed, 25 pieces each stacked inside each other, 500 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

## Size

Length excluding cap	71 mm
Length including cap	73 mm
Diameter of opening	58 mm
Max. volume	100 ml
Sample volume	100 ml

## Product characteristics

Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	12000
Height of inner box	325 mm
Width of inner box	404 mm
Depth of inner box	600 mm
Height of case	325 mm
Width of case	404 mm
Depth of case	600 mm
Volume of case	0.0788 cbm
Weight of case	7.64 kg
EAN inner box	4038917247801
EAN case	8052109529632

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT