



## Product specification

### miniPERM®, Disposable syringe 2 ml, for miniPERM® bioreactor, Luer



## Product information

|                |             |
|----------------|-------------|
| Product number | 94.6077.136 |
|----------------|-------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | miniPERM®, Disposable syringe 2 ml,<br>for miniPERM® bioreactor, Luer   |
| Product number      | 94.6077.136   |
| Product description | miniPERM®, Disposable syringe 2 ml,<br>for miniPERM® bioreactor, Luer, not<br>autoclavable, sterile, 1 piece(s)/blister |

## Size

|             |      |
|-------------|------|
| Max. volume | 2 ml |
|-------------|------|

## Product characteristics

|            |      |
|------------|------|
| Connection | Luer |
|------------|------|

## Purity & certification

|                  |                     |
|------------------|---------------------|
| Product category | Medical device (MD) |
| CE certified     | CE supplier         |
| Purity standard  | sterile             |
| Sterilisation    | Gas sterilisation   |
| Pyrogen-free     | yes                 |

Batched

yes

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 100           |
| Piece(s) / case      | 100           |
| Piece(s) / pallet    | 32400         |
| Height of case       | 100 mm        |
| Width of case        | 191 mm        |
| Depth of case        | 268 mm        |
| Volume of case       | 0.0051 cbm    |
| Weight of case       | 0.42 kg       |
| EAN inner box        | 4038917255929 |
| EAN case             | 4038917255912 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT