



Product specification

**miniPERM®, Disposable syringe 50 ml, for miniPERM®
bioreactor, Luer lock**



Product information

Product number	94.6077.137
----------------	-------------

Product descriptions

Name	miniPERM®, Disposable syringe 50 ml, for miniPERM® bioreactor, Luer lock
------	--

Product number	94.6077.137
----------------	-------------

Product description	miniPERM®, Disposable syringe 50 ml, for miniPERM® bioreactor, Luer lock, not autoclavable, sterile, 1 piece(s)/blister
---------------------	---

Size

Max. volume	50 ml
-------------	-------

Product characteristics

Connection	Luer lock
------------	-----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE supplier
--------------	-------------

Purity standard	sterile
-----------------	---------

Sterilisation	Electron irradiation
---------------	----------------------

Batched	yes
---------	-----

Packaging

Piece(s) / case	60
Piece(s) / pallet	3840
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.53 kg
EAN inner box	4038917255950
EAN case	4038917255943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT