



Product specification

Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 10.8 mm



Product information

Product number	95.64.249
----------------	-----------

Product descriptions

Name	Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 10.8 mm
Product number	95.64.249
Product description	Rack, lid enclosed, material: styrofoam, white, format: 20 x 5, (LxWxH): 340 x 100 x 65 mm, for 100 collection tubes, suitable for tubes Ø 10.8 mm, 1 piece(s)/blister

Material & colours

Color of product	white
Closure material	Styrofoam
Color of cap	white

Size

Width of product	100 mm
Length of product	340 mm
Height of product	65 mm

Storage capacity	100
------------------	-----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	30
Piece(s) / pallet	480
Height of inner box	68 mm
Width of inner box	102 mm
Depth of inner box	343 mm
Height of case	380 mm
Width of case	428 mm
Depth of case	498 mm
Volume of case	0.081 cbm
Weight of case	2.3 kg
EAN inner box	4038917095570
EAN case	4038917007757

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT