



Product specification

PCR film, free from DNase/RNase, transparent



Product information

Product number	95.1992
----------------	---------

Product descriptions

Name	PCR film, free from DNase/RNase, transparent
------	--

Product number	95.1992
----------------	---------

Product description	PCR film, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), free from DNase/RNase, transparent, suitable for qPCR, self-adhesive, application area: -80 - 120 °C, 100 piece(s)/box
---------------------	---

Material & colours

Color of product	transparent
------------------	-------------

Product characteristics

Light protection	no
------------------	----

Closure type	Film
--------------	------

Pierceable	no
------------	----

Suitable for qPCR	yes
-------------------	-----

Minimum temperature	-80 °C
---------------------	--------

Maximum temperature	120 °C
---------------------	--------

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNase/RNase
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	90000
Height of inner box	170 mm
Width of inner box	46 mm
Depth of inner box	90 mm
Height of case	380 mm
Width of case	380 mm
Depth of case	380 mm
Volume of case	0.0549 cbm
Weight of case	20 kg
EAN inner box	4038917940078
EAN case	4038917940085

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT