



Product specification

Storage box, slip-on lid, cardboard, format: 12 x 12, for 144 collection tubes



Product information

Product number	95.64.924
----------------	-----------

Product descriptions

Name	Storage box, Slip-on lid, Cardboard, Format: 12 x 12, for 144 Collection tubes
Product number	95.64.924
Product description	Storage box, slip-on lid, material: cardboard, white, format: 12 x 12, for 144 collection tubes, suitable for tubes up to max. 9.5 mm Ø, tube height between 19-26 mm, 1 piece(s)/bag

Material & colours

Color of product	white
Closure material	Cardboard
Color of cap	white

Size

Width of product	135 mm
Length of product	135 mm
Height of product	28 mm

Product characteristics

Storage capacity	144
Closure type	Slip-on lid

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / case	216
Piece(s) / pallet	2592
Height of inner box	27 mm
Width of inner box	135 mm
Depth of inner box	135 mm
Height of case	521 mm
Width of case	521 mm
Depth of case	610 mm
Volume of case	0.1653 cbm
Weight of case	15.6 kg
EAN inner box	4038917095723
EAN case	4038917067461

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT