




## Product specification

**light blue, suitable for S-Monovette®, tubes Ø 13-16 mm**



### Product information

Product number	65.647.320
----------------	------------

### Product descriptions

Name	light blue, suitable for S-Monovette®, tubes Ø 13-16 mm
Product number	65.647.320
Product description	with filter, light blue, suitable for S-Monovette®, tubes Ø 13-16 mm, 1,000 piece(s)/bag

### Material & colours

Color of product	light blue
------------------	------------

### Product characteristics

Filter	yes
--------	-----

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Piece(s) / inner box	1000
Piece(s) / case	5000

Piece(s) / pallet	210000
Height of case	388 mm
Width of case	250 mm
Depth of case	378 mm
Volume of case	0.0367 cbm
Weight of case	5.37 kg
EAN inner box	4038917988278
EAN case	4038917988285

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT