



Product specification

Filtropur V100, Vacuum filtration unit, 1,000 ml, PES, 0.22 µm



Product information

Product number	83.3942.501
----------------	-------------

Product descriptions

Name	Filtropur V100, Vacuum filtration unit, 1,000 ml, PES, 0.22 µm
Product number	83.3942.501
Product description	Filtropur V100, Vacuum filtration unit, 1,000 ml, screw cap, membrane: PES, Ø membraner: 91 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	High Density Polyethylene (HD-PE)

Size

Height of product	31 cm
Sample volume	1000 ml
Filtration volume	1000 ml

Product characteristics

Filter	yes
Area of the filter	58.93 cm ²
Graduation	yes
Pore size	0.22 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for sterile filtration

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	12
Piece(s) / pallet	192
Height of case	350 mm
Width of case	425 mm
Depth of case	545 mm
Volume of case	0.0811 cbm
Weight of case	5.42 kg
EAN inner box	4038917969123
EAN case	4038917969130

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT