



Product specification

Polystyrene container, suitable as outer packaging for cool transport



Product information

Product number	95.1011
----------------	---------

Product descriptions

Name	Polystyrene container, suitable as Outer packaging for cool transport
Product number	95.1011
Product description	Polystyrene container, suitable as outer packaging for cool transport

Material & colours

Color of product	white
------------------	-------

Size

Width of product	71 mm
Length of product	140 mm
Height of product	168 mm

Product characteristics

Design	Outer packaging for cool transport
Suction insert	no

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	36
Piece(s) / pallet	720
Height of case	350 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0791 cbm
Weight of case	2.25 kg
EAN inner box	4038917084987
EAN case	4038917362399

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT