



Product specification

Blood gas capillaries, calcium-balanced lithium heparin, 100 µl



Product information

Product number	19.930.100
----------------	------------

Product descriptions

Name	Blood gas capillaries, calcium-balanced lithium heparin, 100 µl
Product number	19.930.100
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 100 µl, (LxØ): 100 x 2.05 mm, 250 piece(s)/container

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	2.05 mm
Length of product	100 mm
Sample volume	100 µl

Product characteristics

Preparation concentration	80 IU/ml blood
Type of collection	capillary

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Type of inner packaging	Case
Piece(s) / inner box	250
Piece(s) / case	16000
Piece(s) / pallet	960000
Height of case	175 mm
Width of case	358 mm
Depth of case	398 mm
Volume of case	0.0249 cbm
Weight of case	6.336 kg
EAN inner box	4038917234252
EAN case	4038917234245

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT