



Product specification

Blood gas capillaries, calcium-balanced lithium heparin, 125 µl



Product information

Product number	19.930.125
----------------	------------

Product descriptions

Name	Blood gas capillaries, calcium-balanced lithium heparin, 125 µl
Product number	19.930.125
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	2.05 mm
Length of product	110 mm
Sample volume	125 µl

Product characteristics

Preparation concentration	80 IU/ml blood
Type of collection	capillary

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Type of inner packaging	Case
Piece(s) / inner box	250
Piece(s) / case	16000
Piece(s) / pallet	1280000
Height of case	180 mm
Width of case	288 mm
Depth of case	378 mm
Volume of case	0.0196 cbm
Weight of case	6.64 kg
EAN inner box	4038917234283
EAN case	4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT