



Product specification

Blood gas capillaries, calcium-balanced lithium heparin, 140 µl



Product information

Product number	19.930.140
----------------	------------

Product descriptions

Name	Blood gas capillaries, calcium-balanced lithium heparin, 140 µl
Product number	19.930.140
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 140 µl, (LxØ): 75 x 2.3 mm, 200 piece(s)/container

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	2.3 mm
Length of product	75 mm
Sample volume	140 µl

Product characteristics

Preparation concentration	80 IU/ml blood
Type of collection	capillary

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Type of inner packaging	Case
Piece(s) / inner box	200
Piece(s) / case	12800
Piece(s) / pallet	921600
Height of case	190 mm
Width of case	303 mm
Depth of case	333 mm
Volume of case	0.0192 cbm
Weight of case	4.224 kg
EAN inner box	4038917234313
EAN case	4038917234306

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT