



Product specification

**UriSet 24, without preparation, 3 l, with Urine-Monovette®
10 ml, brown, with light protection, grad**



Product information

Product number	77.578.010
----------------	------------

Product descriptions

Name	UriSet 24, without preparation, 3 l, with Urine-Monovette® 10 ml, brown, with light protection, graduated
------	---

Product number	77.578.010
----------------	------------

Product description	UriSet 24, Complete set for optimum 24-hour urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, without preparation, brown, with light protection, graduated, material: PE, cap: yellow, with Urine-Monovette® 10 ml, 30 piece(s)/case
---------------------	--

Material & colours

Color of product	brown
------------------	-------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	yellow
--------------	--------

Size

Length excluding cap	243 mm
----------------------	--------

Length including cap	246 mm
----------------------	--------

Width of product	130 mm
------------------	--------

Length of product	125 mm
Diameter of opening	83.5 mm
Max. volume	3 l
Sample volume	3000 ml

Product characteristics

Design	with Urine-Monovette® 10 ml
Graduation	yes
Light protection	yes
Viewing window	yes
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	270
Height of case	515 mm
Width of case	403 mm
Depth of case	643 mm
Volume of case	0.1335 cbm
Weight of case	6.34 kg
EAN inner box	4038917113861
EAN case	4038917080019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



SARSTEDT

This document was prepared by EDP support and is valid without signature.