




Product specification

Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 13 mm



Product information

Product number	95.64.213
----------------	-----------

Product descriptions

Name	Rack, styrofoam, format: 20 x 5, suitable for tubes Ø 13 mm
Product number	95.64.213
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 336 x 89 x 49 mm, for 100 collection tubes, suitable for tubes Ø 13 mm, 1 piece(s)/blister

Material & colours

Color of product	white
Closure material	Styrofoam
Color of cap	white

Size

Width of product	89 mm
Length of product	336 mm
Height of product	49 mm

Product characteristics

Storage capacity	100
------------------	-----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	35
Piece(s) / pallet	700
Height of case	345 mm
Width of case	405 mm
Depth of case	455 mm
Volume of case	0.0636 cbm
Weight of case	1.75 kg
EAN inner box	4038917095563
EAN case	4038917030809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT