



Product specification

Rack, styrofoam, format: 5 x 5



Product information

Product number	95.64.257
----------------	-----------

Product descriptions

Name	Rack, styrofoam, format: 5 x 5
Product number	95.64.257
Product description	Rack, material: styrofoam, white, format: 5 x 5, (LxWxH): 204 x 167 x 32 mm, for 25 collection tubes

Material & colours

Color of product	white
------------------	-------

Size

Width of product	167 mm
Length of product	204 mm
Height of product	32 mm

Product characteristics

Storage capacity	25
------------------	----

Packaging

Piece(s) / inner box	24
Piece(s) / case	24
Piece(s) / pallet	864

Height of case	258 mm
Width of case	464 mm
Depth of case	413 mm
Volume of case	0.0492 cbm
Weight of case	0.5 kg
EAN inner box	4038917095624
EAN case	4038917007801

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT