



Product specification

Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm



Product information

Product number	95.937
----------------	--------

Product descriptions

Name	Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm
Product number	95.937
Product description	Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm, blood gas analysis, 10 piece(s)/case

Material & colours

Color of product	red
------------------	-----

Size

Diameter	19 mm
Length of product	25 mm
Height of product	25 mm

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Type of inner packaging	Case
Piece(s) / inner box	10
Piece(s) / case	160
Piece(s) / pallet	51840
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	1.2901 kg
EAN inner box	4038917234641
EAN case	4038917234634

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT