



Product specification

SafeSeal reaction tube, 5 ml, PP, PCR Performance

Tested, Low DNA-binding



Product information

Product number	72.701.700
----------------	------------

Product descriptions

Name	SafeSeal reaction tube, 5 ml, PP, PCR Performance Tested, Low DNA-binding
------	---

Product number	72.701.700
----------------	------------

Product description	SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low DNA-binding, 25 piece(s)/bag
---------------------	---

Material & colours

Color of product	transparent
------------------	-------------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	natural
--------------	---------

Size

Diameter	17 mm
----------	-------

Length of product	60 mm
-------------------	-------

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	25000 x g

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	200
Piece(s) / pallet	15600
Height of inner box	97 mm
Width of inner box	151 mm
Depth of inner box	189 mm
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	1.1238 kg
EAN inner box	4038917905992
EAN case	4038917906005

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT