



## Product specification

### Screw cap micro tube, 0.5 ml, sterile



## Product information

Product number	72.730.311
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 0.5 ml, sterile
Product number	72.730.311
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	2500
Piece(s) / pallet	130000
Height of case	135 mm
Width of case	545 mm
Depth of case	375 mm
Volume of case	0.0276 cbm
Weight of case	3.58 kg
EAN inner box	4038917994576
EAN case	4038917994583

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT