



Product specification

Deep Well plate, Low protein-binding, 2.2 ml, PCR

Performance Tested, PP



Product information

Product number	82.1972.600
----------------	-------------

Product descriptions

Name	Deep Well plate, Low protein-binding, 2.2 ml, PCR Performance Tested, PP
Product number	82.1972.600
Product description	Deep Well plate, Low protein-binding, 96 well, 2.2 ml, compatible with KingFisher™ Flex/Duo Prime/Presto/Apex, Bio Sprint 96, Chemagic™ Prime™, STARlet, MultiMACS M96 Separator, conical wedge, without label, PCR Performance Tested, material: PP, Edged cavities, 4 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	7.52 mm
Width of product	85 mm
Length of product	127.76 mm
Height of product	41.85 mm
Max. volume	2.2 ml

Sample volume

2.2 ml

Product characteristics

Number of wells	96
Base shape	Conical wedge
Label/Print	without label
Closure type	Film
Centrifugation max (RCF)	4500 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	2304
Height of case	210 mm
Width of case	268 mm
Depth of case	283 mm
Volume of case	0.0159 cbm
Weight of case	2.4284 kg
EAN inner box	4038917957892
EAN case	4038917957908

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT