

# Product specification

## Storage box, slip-on lid, cardboard, format: 7 x 7, for 49 collection tubes



### Product information

Product number	95.64.950
----------------	-----------

### Product descriptions

Name	Storage box, Slip-on lid, Cardboard, Format: 7 x 7, for 49 Collection tubes
------	---

Product number	95.64.950
----------------	-----------

Product description	Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes Ø 16-17 mm
---------------------	---

### Material & colours

Color of product	white
------------------	-------

Closure material	Cardboard
------------------	-----------

Color of cap	white
--------------	-------

### Size

Width of product	135 mm
------------------	--------

Length of product	135 mm
-------------------	--------

Height of product	87 mm
-------------------	-------

### Product characteristics

Storage capacity	49
------------------	----

Closure type

Slip-on lid

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / case	60
Piece(s) / pallet	960
Volume of case	0.1134 cbm
Weight of case	8.27 kg
EAN inner box	4038917319607
EAN case	4038917319591

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

