



Product specification

Storage box, slip-on lid, cardboard



Product information

Product number	95.64.998
----------------	-----------

Product descriptions

Name	Storage box, Slip-on lid, Cardboard
Product number	95.64.998
Product description	Storage box, slip-on lid, material: cardboard, white

Material & colours

Color of product	white
Closure material	Cardboard
Color of cap	white

Size

Width of product	135 mm
Length of product	135 mm
Height of product	45 mm

Product characteristics

Label/Print	With print
Closure type	Slip-on lid

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched

yes

Packaging

Piece(s) / inner box	10
Piece(s) / case	120
Piece(s) / pallet	2160
Volume of case	0.1134 cbm
Weight of case	7.86 kg
EAN inner box	4038917246927
EAN case	4038917246910

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT