



Product specification

Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm



Product information

Product number	95.936
----------------	--------

Product descriptions

Name	Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm
Product number	95.936
Product description	Mixing wire for blood gas capillaries, (LxØ): 9 x 1 mm, accessories for the Blood Gas Monovette®, 250 piece(s)/container

Size

Diameter	1 mm
Length of product	9 mm

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Type of inner packaging	Rack
-------------------------	------

Piece(s) / inner box	250
Piece(s) / case	41250
Piece(s) / pallet	412500
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	2.4338 kg
EAN inner box	4038917234610
EAN case	4038917234603

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT