



Product specification

Post transport packaging, 192 x 146 x 75 mm, for cold transport container



Product information

Product number	95.903
----------------	--------

Product descriptions

Name	Post transport packaging, 192 x 146 x 75 mm, for Cold transport container
Product number	95.903
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, (LxWxH): 192 x 146 x 75 mm, for cold transport container, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, 50 piece(s)/case

Material & colours

Product material	Cardboard
Color of product	white

Size

Width of product	146 mm
Length of product	192 mm
Height of product	75 mm

Product characteristics

Label/Print	With print
Color of print/label	black
Example of use	External packaging according to P 650/ ADR

Purity & certification

Satisfies the requirement	P 650, ADR, IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	1650
Height of case	115 mm
Width of case	513 mm
Depth of case	548 mm
Volume of case	0.0323 cbm
Weight of case	4.54 kg
EAN inner box	4038917084581
EAN case	4038917079594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT