



## Product specification

Post transport packaging, 192 x 146 x 77 mm, for cold transport container



## Product information

Product number	95.903
----------------	--------

## Product descriptions

Name	Post transport packaging, 192 x 146 x 77 mm, for cold transport container
Product number	95.903
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, (LxWxH): 192 x 146 x 77 mm, for cold transport container, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, 50 piece(s)/case

## Material & colours

Color of product	white
------------------	-------

## Size

Width of product	146 mm
Length of product	192 mm
Height of product	77 mm

## Product characteristics

Label/Print	With print
-------------	------------

Color of print/label	black
----------------------	-------

Example of use	External packaging according to P 650/ ADR
----------------	--

## Purity & certification

Satisfies the requirement	P 650, ADR, IATA
---------------------------	------------------

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Batched	yes
---------	-----

## Packaging

Piece(s) / inner box	50
----------------------	----

Piece(s) / case	50
-----------------	----

Piece(s) / pallet	1650
-------------------	------

Height of case	115 mm
----------------	--------

Width of case	513 mm
---------------	--------

Depth of case	548 mm
---------------	--------

Volume of case	0.0323 cbm
----------------	------------

Weight of case	4.54 kg
----------------	---------

EAN inner box	4038917084581
---------------	---------------

EAN case	4038917079594
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT