



Product specification

**PCR plate half skirt, 96 well, white, Low Profile, 100 µl,
PCR Performance Tested, PP**



Product information

| | |
|----------------|-------------|
| Product number | 72.1982.252 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|---|
| Name | PCR plate half skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP |
| Product number | 72.1982.252 |
| Product description | PCR plate half skirt, for LightCycler® 480, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 10 piece(s)/minigrip bag |

Material & colours

| | |
|------------------|-------------|
| Color of product | white |
| Profile | Low Profile |

Size

| | |
|-------------------|----------|
| Diameter | 5.38 mm |
| Width of product | 85.3 mm |
| Length of product | 127.8 mm |
| Height of product | 15.55 mm |
| Max. volume | 100 µl |
| Sample volume | 0.1 ml |

Product characteristics

| | |
|--|------------|
| Number of wells | 96 |
| Construction of the edges of the plate | half skirt |
| Base shape | conical |
| Closure type | Film |
| Autoclavable (max. °C) | 121 °C |
| Suitable for qPCR | yes |
| Centrifugation max (RCF) | 3000 x g |
| Centrifugation temperature | 22 °C |
| Centrifugation time | 5 Min |

Purity & certification

| | |
|------------------|------------------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 50 |
| Piece(s) / case | 100 |
| Piece(s) / pallet | 7200 |
| Height of inner box | 140 mm |
| Width of inner box | 170 mm |
| Depth of inner box | 310 mm |
| Height of case | 190 mm |
| Width of case | 303 mm |
| Depth of case | 333 mm |
| Volume of case | 0.0192 cbm |
| Weight of case | 3.213 kg |
| EAN inner box | 4038917932899 |
| EAN case | 4038917932905 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT